



Fetal monitoring device

FC1400

Operation Manual



Ver. 2.42

2023.04.10



www.ebionet.com



REVISION HISTORY

[illegible]

Warranty

- This product is manufactured through our strict quality control and inspection process. Compensation standards for product repair and exchange follow the "Regulations of Compensation for Consumer's Damage" announced by the Fair-Trade Commission.
- Warranty period of this product is regulated to be 1 year while the warranty period of accessories is six months.
- If a malfunction occurs under normal use, our service center will repair it free of charge during the warranty period.
- If a problem occurs with the product during the warranty period, please notify us of the model name, serial number, date of purchase, and malfunction details.

CAUTION
Federal law restricts this device to sale by or on the order of a physician

NOTE
The product does not have shelf life. Its expected use life is 6.5 years. After 6.5 years, though the product still works normally, it is recommended to have it checked by Bionet.

Contact Bionet

If you have any questions or comments relating to our products or purchasing, please contact the telephone numbers or E-mail below. You can talk to our sales people.

Bionet always welcomes your enquiries. Please contact us.

**Headquarters &
International****Sales & service****Bionet Co., Ltd.:**

5F, 61 Digital-ro 31 gil, Guro-gu, SEOUL 08375,
REPUBLIC OF KOREA

Tel : +82-2-6292-6410 / Fax : +82-2-6499-7789

E-mail : service@ebionet.com

Website: www.ebionet.com

**U.S.A
sales & service
representative****Bionet America, Inc.**

2691 Dow Ave. Suite B, Tustin, CA 92780

Toll Free : 1-877-924-6638 / Fax : 1-714-734-1761

E-mail : support@bionetus.com

Website: www.bionetus.com

**European
sales & service
representative****Bionet Europe GmbH :**

2Li Bessemerstr. 51, D-12103 Berlin, Germany

Tel : +49-30-240-374-52

E-mail : bionetEU@ebionet.com

Website: www.bioneteu.com

**European
Representative****CMC Medical Devices & Drugs S.L. :**

C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain

Tel +34-951-214-054 / Fax +34-952-330-100

E-mail: info@cmcmedicaldevices.com

Website: www.cmcmedicaldevices.com

※ In the event of a malfunction or failure, contact Service Dept. Of Bionet Co., Ltd.
along with the model name, serial number, date of purchase and explanation of failure.

Paid Services

A fee will be charged for all services except for breakdowns, so be sure to read this Operation manual below before putting in a request.

- Usage description and simple inspection without disassembly - In case of reinstallation due to poor installation by a distributor	Free the 1st time Charged starting the 2nd time
- Inadequate installation or loosening due to physical product movement, relocation, etc. - When re-installing after the first installation requested by the customer at the time of purchase - When reinstallation is required due to inexperienced installation by the customer - When a service is requested due to the input of foreign substances or improper cleaning	Charged starting the 1st time

1. Equipment cleaning, adjustment, and usage description are not product breakdowns.

(Unfeasible repairs are subject to separate standards.)

2. Breakdowns caused by consumer negligence

Breakdowns and damage due to careless handling by the customer or incorrect repair are caused by:

- Using incompatible electric capacity.
- Mishaps such as dropping the product.
- Using the third party replacements or options not specified by our company.
- Non-Bionet technicians or agency technicians in the process of repair.

3. Other cases

- Breakdowns by natural disasters (fire, salt damage, flood damage, earthquake, etc.)
- When a consumable part has reached the end of its life (accessories)

Table of Contents

CHAPTER 1. BASICS	9
DEFINITION OF WARNING, PROHIBITION, MANDATORY ACTION, AND NOTE	9
CAUTIONS.....	11
BIOCOMPATIBILITY	16
MAINTENANCE, CLEANING AND EQUIPMENT CONNECTION	17
SAFETY SYMBOLS	20
CHAPTER 2. INSTALLATION	22
1) PRODUCT OVERVIEW	22
2) PRODUCT CHARACTERISTICS.....	22
3) COMPOSITION OF PRODUCTS.....	23
4) INSTALLING SYSTEM.....	30
CHAPTER 3. BASIC OPERATIONS	33
1) GETTING STARTED	33
2) SHUTTING DOWN	33
3) GRAPHIC DISPLAY	33
4) CONTROL PANEL.....	35
5) MENU.....	37
6) CONNECTING POWER	38
7) DEFINITION OF GROUPS	41
CHAPTER 4. PATIENT MANAGEMENT	42
1) REGISTERING A NEW PATIENT	42
2) EDITING A PATIENT.....	44
CHAPTER 5. MEASURING FHR USING US	45
1) MEASURING FHR.....	45
2) US SCREEN	48
3) SETUP	48
CHAPTER 6. MEASURING THE UTERINE CONTRACTION EXTERNALLY	51
1) MEASURING UTERINE CONTRACTION	51
2) UC SCREEN	53
3) SETUP	53
CHAPTER 7. MEASURING FETAL MOVEMENT.....	55
CHAPTER 8. CONNECTING THE STIMULATOR.....	57

CHAPTER 9. CLINICAL MARKS AND NOTES.....	58
CHAPTER 10. TRACE	60
1) TRACE AREA.....	60
2) TRACE.....	61
3) TRACE MENU	62
CHAPTER 11. PRINTING	66
1) REAL-TIME PRINTING	67
2) PRINTING THE TRACE	68
3) LOADING THE PAPER.....	68
4) REGISTERING HOSPITAL NAME	68
CHAPTER 12. ALARM AND PRESET.....	69
1) PATIENT CONDITION ALARMS.....	69
2) PRODUCT CONDITION ALARMS.....	69
3) VISUAL ALARMS	70
4) ALARM LED	71
5) AUDIBLE SILENCE AND PAUSE.....	71
6) ALARM LATCHING.....	72
7) ALARM HISTORY	72
8) PRESET	72
9) ADJUSTING ALARM VOLUME	73
10) SETTING ALL ALARMS ON/OFF	74
11) SETTING RANGES OF ALL ALARMS	75
11) SETTING LEVELS OF ALL ALARMS.....	75
13) DEFAULT SETTING.....	76
14) SIGNAL LOSS ALARM.....	76
CHAPTER 13. NETWORK	77
1) SETTING NETWORK.....	77
CHAPTER 14. GENERAL SETTINGS	80
1) CHANGING THE DATE	80
2) CHANGING THE TIME	80
3) CHANGING THE LANGUAGE	80
4) CHECKING THE SYSTEM VERSION	81
5) CHANGING THE TOUCH TONE.....	81
6) DEMO OPERATION	81
7) EDITING NOTE	82

8) MARKER SOUND	82
9) FACTORY MODE.....	82
10) CHANGING ADMIN PW	82
11) VOLUME RANGE.....	82
12) SCREEN OUTPUT MODE	82
13) PROTOCOL VERSION	83
CHAPTER 15. NON STRESS TEST (NST).....	84
1) MEASURING NST	84
2) SETTING NST	84
CHAPTER 16. CTG	85
1) SETUP	85
2) PRINTING CTG	85
3) CTG MEASUREMENT RESULTS.....	86
4) GLOSSARY OF CTG TERMS	87
CHAPTER 17. MESSAGE LIST	90
CHAPTER 18. SIMPLE TROUBLESHOOTING	91
1) TROUBLESHOOTING AND SOLUTIONS.....	91
2) PERFORMING PERIODIC INSPECTIONS	91
3) CYBER SECURITY ISSUES	91
CHAPTER 19. PRODUCT SPECIFICATIONS	92
APPENDIX A.	95
APPENDIX B. MAINTENANCE, CARE, AND SERVICE	101
APPENDIX C. ULTRASOUND POWER	104

Specifications or functions described in this operation manual are subject to change without notice for product improvement.

Chapter 1. Basics

Definition of Warning, Prohibition, Mandatory Action, and Note

- The following terms are used throughout this manual to emphasize important and critical information. You must read these statements to help ensure safety and to prevent product damage.
- The manufacturer or the product distributor is not liable for any loss or damage to the product caused by incorrect use or negligence in product maintenance.

WARNING	
	Failure to follow this message may cause severe injuries, casualty or physical damage to patients.










CAUTION	
	Failure to follow this message may cause in non-life-threatening injury or damage to the equipment.

MANDATORY ACTION	
	Mandatory Action messages are directions you must follow for safe operation and maintenance of the equipment.

NOTE	
Note messages indicate some important information and tips, which are not dangerous, about installation, operations and maintenance.	

General Precautions on Environment

DO NOT store or operate the equipment in the places listed below.

	A place exposed to moisture (DO NOT touch the equipment with wet hands.)		A place under direct sunlight
	A place in areas with highly fluctuating temperatures.		A place in the vicinity of Electric heater
	A place with excessive humidity rise or poor ventilation		A place with sources that cause excessive shock or vibration
	A place exposed to chemicals or at risk of gas leakage		Avoid the invasion of small objects/ particles such as dust, and especially avoid metallic material.
	DO NOT disjoint or disassemble the equipment. (Bionet is not liable for broken products caused by attempted disassembly.)		DO NOT connect power until the product is completely installed. It may cause damage to the product.

CAUTIONS

Before Installation

Compatibility is critical to safe and effective use of this equipment. Please contact your local sales or service representative prior to installation to verify equipment compatibility.


Defibrillator Precautions

The patient signal inputs labeled with the CF and BF symbols are protected against damage resulting from defibrillation voltages. To ensure proper defibrillator protection, use only the recommended cables and lead wires. Proper placement of defibrillator paddles in relation to the electrodes is required to ensure successful defibrillation.

Disposables

Disposable devices are intended for single use only. They should not be reused as performance could degrade or contamination could occur.

Existing Device Disposal



1. Products bearing this symbol (X-marked wheeled bins) are subject to European Directive 2002/96/EC.
2. All electrical and electronic products must be disposed of separately from municipal waste at the collection facility designated by government or local authorities.
3. Proper disposal of old devices helps prevent potential adverse consequences against environmental and human health.
4. For more information on the disposal of existing devices, contact City Hall, Waste Disposal Service Center, or the store where you purchased the product.

WARNING
This equipment contains a chemical known to the State of California to cause cancer, birth defects, or other reproductive harm.

Electrocute Precautions

To prevent skin burns, apply electrocute electrodes as far as possible from all other electrodes. A distance of at 15 cm/6 in. is recommended.

EMC

Magnetic and electrical fields are capable of interfering with the proper performance of this equipment. For this reason, make sure that all external devices operated in its vicinity comply with the relevant EMC requirements. X-ray equipment or MRI devices are possible source of interference as they may emit higher levels of electromagnetic radiation. Also, keep cellular phones to other telecommunication equipment away from it.

Instruction for Use

To continue using this equipment safely, it is necessary that you follow the instructions. However, instructions listed in this manual in no way supersede established medical practices concerning patient care.

Loss of Data

Should this equipment at any time temporarily lose patient data, the potential exists that active monitoring is not being done.

Close patient observation or alternate monitoring devices should be used until the monitoring function is restored. If this equipment does not automatically resume operation within 60 seconds, power cycle it using the power on/off switch. Once monitoring is restored, you should verify correct the monitoring state and alarm function.

Maintenance

Regular technical inspections should be carried out annually to meet any requirements specific to your country.

MPSO

The use of a multiple portable socket outlet (MPSO) for this equipment result in current leakage in it. DO NOT use MPSO as much as possible.

Negligence

Bionet Co., Ltd. does not assume responsibility for damage to this equipment caused by improper or faulty power or incorrect placement.

NOTES

Power Requirements

This equipment uses DC adapter of 100-240VAC / 18VDC 2.8A. Be sure to use the adapter provided by Bionet.

Restricted Sale

U.S.A federal law restricts this equipment to sale by hospital or on the order of a physician.

Supervised Use

This equipment is intended for use under the direct supervision of a licensed health care practitioner.

Installation Requirements

Set up this equipment in a location which affords sufficient ventilation.

(The ambient conditions specified in the technical specifications must be ensured at all times.)

Put this equipment in a location where you can easily see the screen and access the operating controls.

Safety Instructions for Electricity

Please note the following precautions before using the product.

- Is the power supply cord proper? (100 - 240V AC)
- Is every cord connected properly to the product?
- Is the product fully grounded? (Otherwise, noise may occur.)
- Is there any damage - which may affect the monitoring or safety of the patients — to any of the product or accessories?

NOTE

You should not place this equipment in the vicinity of electric generator, X-ray, broadcasting apparatus or portable wires to eliminate electric noise during operation. Otherwise, it may cause incorrect results. Isolated power line and stable grounding are important requirements for use of this equipment. Using same power source with other electric equipment may cause incorrect results.

NOTE

FC1400 is classified as listed below:

- As for protection against electric shock, it is a Class 2 type and has type-BF degree.
- Do not use it in the vicinity of flammable anesthetics and solvents.
- It fulfills the level A according to IEC/EN 60601-1-2 (Electromagnetic Compatibility Requirements).

NOTE

Doctors and patients in hospitals are exposed to the risk of uncontrollable currents. These currents are caused by a potential difference between this equipment and conductive objects, which may come into it. To solve this problem, make sure that the auxiliary equipment connected to this equipment satisfies EN60601-1-1.

WARNING



DO NOT make contact with patient during defibrillation. Otherwise, severe injury or death could result. To avoid the risk of serious electrical burn, shock, or other injury during defibrillation, all persons must keep clear of the bed and must not touch the patient, or any devices connected to the patient.

NOTE

If a multiple socket outlet is used for AC power input, it must be a grounding-type.

NOTE

Avoid contact with the connector pin of this equipment and the patient at the same time while it is in use.

Do not connect or remove the power cable with wet hands.

NOTE

In case the medical equipment does not operate normally, or if it has been damaged, do not use it on any patient; Contact the medical equipment technician in your hospital or the equipment supplier.

WARNING

Devices for fetal heartbeat sound measurements using ultrasonic (US) waves or uterine contraction measurements (TOCO) performed from outside of the uterus are not designed for use during any electric operation, defibrillation, or defibrillator discharge.

WARNING

This equipment has not been designed for use with other types of monitoring equipment apart from those devices permitted to use together with it in this manual.

WARNING

Do not touch signal input, signal output or other connectors, and this equipment simultaneously.
Do not connect or remove the power cable with wet hands.

WARNING

This equipment is not designed for use with devices that generate high electromagnetic fields, such as high frequency surgery machines. Therefore, measurements may be affected in the presence of strong electromagnetic sources such as electrosurgery equipment.

Biocompatibility

When used as intended, the parts of this equipment described in this operation manual, including accessories that come in contact with the patient during the intended use, fulfill the biocompatibility requirements of the applicable standards. If you have questions about this matter, please contact Bionet or its representatives.

Maintenance, Cleaning and Equipment Connection

Although FC1400 fetal monitoring device and its accessories can be cleaned in many ways, please use the methods recommended below to avoid unnecessary damage to or contamination of them.

If any dangerous material other than those designated for cleaning has been used, the resulting contaminated or damaged monitor or accessories will not be repaired free of charge even during the warranty period.

Make sure that the monitor, probe, cable, and accessories are free of dirt or dust. Carefully check them after each cleaning or disinfection. If degeneration or damage has been found, do not use the them.

Please take note of the following:

- Be sure not to leave any cleaning/disinfecting chemical residue on the surface of the monitor and accessories. After allowing sufficient time for the chemical to work, wipe off all residues with a cloth damp with water.
- Prevent any fluid from seeping into the monitor, module, or accessories.
- The monitor, module, or accessories should not be soaked in any fluid; protect them from water drops or prevent water from being splashed on them.
- Never use any abrasive material (steel wool or silver polish).
- Never use any bleaching agent.
- Do not use any kind of drying equipment such as heater, oven (including microwave oven), hair dryer, or heating lamp.

NOTE
After cleaning, carefully check the monitor and the probe.

Cleaning Components - Electric cables and lead wires

NOTE
- Do not use acetone or ketone solvents for cleaning. - Do not use an autoclave or steam cleaner. - Do not mix the cleaner with any disinfecting solution as toxic gases may be produced. - Turn off the power before cleaning; do not pour water into the components being cleaned.

Electric cables and lead wires can be wiped off or cleaned with towels wet with lukewarm water,

neutral soap, or isopropyl alcohol. Using ethylene oxide for intensive disinfection (almost complete sterilization) is allowed. Note, however, that it will reduce the lifetime of cables or lead wires. Clean the belt using soap and water while making sure that the water temperature does not exceed 60°C.

NOTE

The decision to sterilize must be made per your institution's requirements with an awareness of the effect on the integrity of the cable or lead wires.

NOTE

FC1400 needs safety inspection once a year. Refer to this operation manual or service manual for inspection items.

After cleaning, carefully inspect the monitor and sensors. Do not use them if they have been damaged or deteriorated.

Clean the exterior of the monitor at least once a month using soft cloth wet with tepid water or alcohol. Do not use lacquer, thinner, ethylene, or any oxidizing agent that may cause damage to it. After verifying that there is no dust or contamination on the cables and accessories, wipe them off with soft cloth wet with 40°C/104°F water. Wipe them at least once a week using clinical alcohol.

NOTE

There is back-up battery on board inside FC1400.
Please dispose of it according to all local regulations.

WARNING

Before changing battery, remove the AC power and check the electrodes of the battery.

If the installation or arrangement of the external grounding wire is doubtful, operate FC1400 with internal power.

If FC1400 is not used for a certain period, remove the backup battery to avoid any safety hazard.

NOTE

Annual Servicing – For continued safety and performance of FC1400, it is recommended that its calibration, accuracy, and electrical safety be verified on an annual basis by a Bionet Service Representative.

Daily Testing – It is essential that FC1400 and its accessories be inspected every day. It is recommended practice to initiate its self-test feature at the beginning of each monitoring session; follow the instructions in Chapter 1 and 2.

Safety Symbols Marked on the Package

Symbols	Contents
	This way up
	Keep dry
	Fragile
	Use no hooks

Safety Symbols

The International Electro Technical Commission (IEC) has established a set of symbols for medical electronic equipment which classify a connection or warn of any potential hazards. The classifications and symbols are shown below

Symbols	Contents
	Power Adapter connection
	External Signal IN/OUT Port
	IEC 60601-1 Type BF applied part
	Power Standby
	USB Port for external signal
	USB Connector
	Conductor provides a connection between the equipment and the potential equalization bus bar of the electrical installation
	Video Out
IPX0, IPX1 and IPX7	Protection against vertically falling water drops (IEC 60529) Water Protection Specification Level 0, Level 1 and Level 7.
	Safety Sign: It indicates that you should read the operation manual. Read this operation manual before starting work or operating the equipment..
	Consult Instructions for Use: This symbol advises the reader to consult the operating instructions for information needed for the proper use of the equipment..

	Waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.
---	--

Chapter 2. Installation

1) Product Overview

FC1400 is a fetal monitoring device used to measure the Fetal Heart Rates (FHR), degree of maternal uterine contraction (UA: Uterine Activity), and fetal movements (FM). It ultrasonically enters the abdomen of the patient and extracts Doppler frequencies from the signals returned after being reflected by the heart of the fetus. The Doppler frequencies vary with the movements of the fetal heart and the changes in the heartbeats of the fetus are output as sounds, which are analyzed to detect the heart rates and movements of the fetus. In addition, it detects the degrees of uterine contraction of the patient using a pressure sensor. It shows the fetal heart rates, maternal uterine activity, and fetal movements on its LCD screen as figures and numeric values, and saves the data in its memory.

2) Product Characteristics

1. FC1400 measures fetal heart rates and fetal movements simultaneously.
2. By using the US probe for twin fetuses, twin fetal FHR and movements can be measured. (optional specification)
3. Measured data can be saved and identified on the LCD. Thus, the condition of the fetus can be efficiently monitored during labor without wasting recording sheets.
4. While reviewing the saved data, any data that should be kept in recording sheets can be printed at high speeds.
5. With ultrasound Doppler probes having high durability even against noise, clear sounds of the heartbeats of the fetus and accurate FHR can be detected.
6. The 7 ultrasonic sensors driven at 1MHz are used to minimize broken FHR waveforms even if the fetus or the patient moves.
7. The fetal movement is automatically detected by the analysis of the Doppler signal, and the intensity and duration are displayed on screen or printed as spike waveforms.
8. Universal communication interfaces including LAN are provided for connection with the central parturition monitoring system.
9. Rechargeable batteries are used; thus, parturition conditions can be continuously monitored even during a power failure. (optional specification)


WARNING



Using non-standard accessories other than the Bionet-provided supplies may cause signal distortion or noise. Be sure to use genuine accessories supplied by Bionet.

3) Composition of Products

The FC1400 fetal monitoring device consists of the following components. Open the packaging box and check that all items below are included. In addition, please check the monitor and the components for any damage).



Basic Composition and Accessories ① Monitor ② US probe (2) ③ UC probe (1) ④ Event Marker (1) ⑤ Fetal Paper (2) ⑥ Power cord (1) ⑦ Adaptor (1) ⑧ Ultrasound gel (1) ⑨ Belt (3) ⑩ FC1400 Stand (1)	Optional specification ⑪ Battery (1)
--	---

WARNING	
	How to replace battery: Please make sure you use the right battery (Lithium-ion battery) as shown here. Otherwise, Bionet is not liable for any damages and/or explosion caused by the wrong battery. Make sure to use the genuine battery recommended by Bionet.  Lithium-ion battery (10.8V / 3250mAh)  Ni-MH Battery (12V / 2600mAh)


WARNING

Do not subject the probe to physical shocks or excessive pressure.


Shock and excessive pressure on the probe may damage the internal sensor.

Impact or excessive pressure on the center of the probe may cause damage.

Bionet is not liable for probes damaged as a result of these actions.



■ Top View




- ① Graphic display window: Displays the measuring items and condition of FC1400.
- ② Alarm LED: Displays the state of alarm.
- ③ Control panel: Controls the functions.
- ④ Printer door: Opens for replacing recording paper.

WARNING




To avoid high-voltage electric shock, DO NOT disassemble the monitor. If any of the probes has been damaged, stop using them to prevent electric shock and contact Bionet or its representative for the replacement of the probe. Disassembling the monitor should be done only by the service representative authorized by Bionet.

■ Front View



- ① Hand Grip: Use it when lifting and moving the monitor.


■ Rear View



- ① Power adaptor connector: Connect a power adaptor of 18V, 2.8A.
- ② LAN Port: A terminal for LAN communication
- ③ USB1 Port: A terminal for USB communication
- ④ USB2 Port: A terminal for USB communication


Mandatory Action	
	DO NOT remove the USB storage device until it is recognized completely.

■ Right Side



- ① Printer door release lever: Slide it to open the printer door.

■ Left Side



- ① Marker jack connection terminal: Connect the marker jack cable.
- ② Stimulator jack connection terminal: Connect the stimulator cable.
- ③ US probe connection terminal 1: A terminal to connect a US probe
- ④ US probe connection terminal 2: A terminal to connect a US probe to measure twin fetal
- ⑤ TOCO probe connection terminal: A terminal to connect the TOCO probe

Prohibition	
	Avoid contact with both the USB or LAN port area and the patient simultaneously. RISK OF ELECTRICAL SHOCKS – DO NOT OPEN THE MONITOR: Monitor disassembly is only authorized to individuals with Bionet service authorization.

4) Installing System

■ Cautions in installation

Preparations before use

- Be sure to use FC1400 in dry condition at normal temperatures (temperature: 10 ~ 40°C, humidity: 30 ~ 85%).
- Be sure to use an outlet away from any equipment that may generate electric noise (room heating/cooling equipment using a large motor, etc.).
- Installing recording paper: Open the printer cover by pushing back the Printer cover open lever on the right side of the monitor. Load paper with printable side up and close the cover.
- Put the adaptor into the connection site at the back of the monitor to connect AC power. At this time, check if the power cord has been installed properly.


Storage Temperature: -20~+60°C (-29°F~+15°F)

Storage Pressure: 50 (500)~106Kpa (1060mbar)

Storage Humidity: 10~95% (without condensation)

Installing FC1400 stand

- Insert the FC1400 stand to the stand hole on the underside of the monitor.



Operating procedure

- Press the power switch on the front side of the monitor for 2 seconds when power has been turned off; on the other hand, press the power switch on the for 3 seconds when power has been turned on.
- Put some ultrasound gel on the US probe, find heartbeat of the fetus and attach the probe around the abdomen of the patient using a fixing belt.
- Fix the TOCO probe on Fundus which is located 10 cm above the belly button of the maternal abdomen using a fixing belt.
- Determine a usage mode and use FC1400 accordingly.
 - * If a problem occurs with FC1400, press and hold the power key for more than 6 seconds to force shutdown.

How to store and manage FC1400 after use

- Turn off the power switch and unplug the power cord.
- When the operation has been completed, clean the monitor and accessories to prevent malfunction.

Cautions in use

Installation and storage

- Avoid moisture, high temperature, dust in the air, salt, and sulfuric materials including places that are under direct sunlight or not well-ventilated.
- Avoid vibrations or mechanical impacts.
- Avoid places exposed to chemicals or at risk of gas leakage.
- FC1400 should be used with the voltage and frequency indicated.

Before operation

- FC1400 should be properly grounded.
- Connect all cords accurately and safely.
- When using an optional battery, please charge it for at least four hours.
- Any part directly connected to the patient should be double-checked.
- Use FC1400 within the range of operating temperatures.

During operation

- The patient should not come into contact with any metal parts of FC1400, such as the chassis or case. DO NOT touch the patient and FC1400 at the same time.
- If any abnormality has been found, turn off power and unplug the power cord.
- DO NOT use any sharp or pointed object to press the LCD or any touch key.
- DO NOT confuse demo data with patient data.

After use

- Turn off the power switch and unplug the power cord.
- Take off the probes connected to the patient.
- When the operation has been completed, clean the monitor and accessories to prevent malfunction.

Periodic Inspection

- Keep the monitor and the measuring probes clean by wiping them off with soft cloth wet with alcohol at least once a month. DO NOT use lacquer, thinner, ethylene, or oxidizing substances.
- Keep the cables free of dust or dirt. Clean them with cloth wet with tepid water (40°C/104°F) after use and at least once a week, wipe them off with clinical alcohol.
- Inspect FC1400 periodically once a year.
- FC1400 should be serviced only by the service representative authorized by Bionet.

Chapter 3. Basic operations

1) Getting Started


Press the power switch on the monitor for 2 seconds; it turns on and home screen appears.

2) Shutting Down

Hold down the power switch for 3 seconds to turn off the monitor. If it does not turn off when you press the power switch for 3 seconds, press and hold the power switch for 6 seconds or more.

3) Graphic Display


■ Monitoring screen



- ① Date and Time
- ② Patient name or ID
- ③ Bed ID
- ④ USB terminal connection
- ⑤ Print state
- ⑥ Network connection
- ⑦ Alarm volume, size, and setting state
- ⑧ Battery capacity
- ⑨ Fetal heart rate parameters
- ⑩ Uterine contraction parameters
- ⑪ Main Menu

- ⑫ The trends of fetal heart rates
- ⑬ The trends in the real-time degrees of uterine contraction and fetal movements

4) Control Panel



- ① AC power indicator LED: AC power connection state
- ② Battery indicator LED: State of battery charges
- ③ US1 volume up/1 min forward key in Trace mode
Short key: US1 Volume up or forwarding 1 minute in trace mode
Long key: Return US1 volume back to the Pre-Mute stage.
- ④ US2 volume up/1 min backward key in Trace mode
Short key: US2 volume up or backwarding 1 minute in Trace mode
Long key: Return US2 volume back to the Pre-Mute stage.
- ⑤ US1 volume down/6 sec forward key in Trace mode
Short key: US1 volume down or forwarding 6 seconds in Trace mode
Long key: US1 volume mute
- ⑥ US2 volume up/6 sec backward key in Trace mode
Short key: US2 volume down or backwarding 6 seconds in Trace mode
Long key: US2 volume mute
- ⑦ Clinical Mark and Note Bottom
Short key: Used by doctors or nurses to mark the current location
Long key: Used by doctors or nurses to record the statements frequently used
- ⑧ Zero calibration uterine contraction: Zero calibration of TOCO probe
- ⑨ Trace entry:

Short key: Used to view trace data.

Long key: Enable or disable the screen touch lock.

- ⑩ Alarm silence key: used to change alarm state. The alarm circulates: alarm on → alarm silence → pause → alarm off → alarm on.
- ⑪ HOME key: You can swap Graphic screen and Text mode. Pressing it in the sub menu moves you to the main menu.
- ⑫ Print start/stop key and Paper feeding key:
 - Short key: Used to print data on the recording paper in real time or in Trace mode.
 - Long key (Paper feeding key): Used to feed papers while printing is stopped.

5) Menu

■ Selecting menu

To go to a menu, touch the main menu or parameter window on the screen.


■ Entering a letter

To enter a letter, touch the key for that character. To clear all entered characters, touch and hold the [←] key for at least 1 second.




■ Entering value

To enter a value, touch the selection key or the screen. If the previously entered value remains, press [CLEAR] key or use [BS] to clear them.



■ Selecting color

To select the color of a line or a numerical value, touch the selection key or color on the screen.



6) Connecting Power

■ AC power


If the monitor is connected to AC power, the Power LED on the front side of the monitor lights green; if a battery has been installed, the power mode is automatically changed to automatic charging mode, and the battery is charged.



■ Battery power (optional)

If you turn on the monitor with the AC power supply being cut off, it is powered by the battery and the battery indicator lights orange in its front. When it is mounted with the battery and connected to the AC power, the battery is charged in automatic charging mode. When the battery is fully charged, the battery indicator LED lights green.

If the battery is low, alarm occurs with *Battery low* being displayed on the LCD screen a few minutes before the monitor is automatically turned off. In this case, connect it to AC power.



For a new battery

- Charging time: Minimum of 4 hours
- Continuous use time: Maximum 2 hours

■ Changes in the screen in relation to the state of battery connection

The battery state is displayed on the screen as follows depending on the connection status of the battery and AC power, and the condition of the battery;

- No battery (AC power connection only)



- Battery only but with little or no power



: 2 images are toggled.

- Battery only and its state



(Low),



(Medium),



(Full)

- AC and Battery: Level of battery being charged.



(Low),



(Medium),




(Full)

■ Replacing the battery

When replacing the battery, use the same type of battery.

- Type:
 - CMICR18650 F9, Lithium Ion rechargeable Battery Pack (10.8V / 3250mV)
- When to replace: When the monitor is connected to a power source, the battery is charged automatically. The battery can be recharged about 300 times and replace it if you cannot use it for longer than 20 minutes even after full charge. In case the battery is damaged or leaks, replace it immediately. Do not use a damaged battery on the monitor.
- How to: Connect the connector as shown in the following image. (Connector cannot be plugged inversely.)



(Connecting a battery)

WARNING

Pay attention to the polarity when replacing the battery.
Be sure to install and remove the battery with the AC power off.

WARNING

To protect the environment, do not carelessly dispose of wastes or residues such as the components of FC1400 after their lifetime; instead, ask the biomedical engineering department of the hospital to dispose of them in the designated places following appropriate procedures.

If the installation or arrangement of the external grounding wire is doubtful, operate FC1400 with internal power. If it is not used for a certain period, remove the battery to avoid any safety hazard.

■ The effect of lithium-ion technology on the battery

Find out about lithium-ion battery technology [here](#).

The battery is discharged naturally even when not installed in the monitor. Discharging is caused by the current demanded by the lithium-ion battery integrated circuit. Battery is self-discharged due to the nature of lithium-ion cells.

Battery retention loss is greater at higher temperatures. As the battery ages, it may not be fully charged, or not charged at all. As a result, the total charge capacity used for saving and using gradually decreases.

Conditioning Guidelines

Check battery performance by fully charging and completely discharging it every 6 months.

Storage Guidelines

Store the battery between 20°C and 25°C (68°F and 77°F) when it is set aside separately. When the battery is installed in the monitor being connected to AC power, the temperature of the battery increases by 15°C to 20°C (59°F to 68°F) at room temperature, which shortens the life of the battery.

When a battery is installed in the monitor with AC power being connected, normally the battery does not use its battery power. Battery life may be less than 12 months. Store the battery along with the monitor to prevent being lost and stolen, and take out the battery from the monitor when moving it around.

7) Definition of Groups

Users of FC1400 are divided into three types, and the authority varies according to each group.

User passwords are group specific.

- Administrator: Can set preset and note items (default password: **1234**)
- Factory: Factory setting. Authorized users at head office and AS representatives
- General user: Can change the settings for the current patient.

Chapter 4. Patient Management

Learn how to store the basic information for patient identification.

1) Registering a New Patient


Register a new patient with this menu.

When you turn on the monitor, a random ID is created with the current date and time.

To change the user's ID to the current ID or to register the changed patient with a different ID, use **New** function. If the records of the patients measured before are left in the list of the recent 99 patients, use **Search** function.

< Registering a new patient >

- ① Touch [Patient] on the menu. (It is not available in Demo mode.)
- ② Select [New] to register a new patient.
- ③ Touch the ID, first name, last name, and gestational-age respectively to enter the text in the input window, then press [OK] to enter the items.
- ④ Enter the ID, first name, last name, and gestational-age and press [OK].
To register the patient with the default values of date and time, touch the [Default] button.
- ⑤ Information on the currently registered patient appears at the patient information tap on top of the screen.



The screenshot shows a 'New Patient' registration form. It has a blue header bar with the title 'New Patient' and a close button (X). The form contains four input fields: 'ID', 'First Name', 'Last Name', and 'Gestational-age'. The 'Gestational-age' field is split into 'weeks' and 'days' sub-fields. At the bottom, there are two buttons: 'Default' and 'Admit'.

< Registering a patient using *Search* >

- ① Touch [Patient] on the menu.
- ② Select [Search] to register a recently measured patient.
- ③ Referring to the ID, first name, last name, and gestational-age on the patient list, touch and select the patient you want to register.
- ④ If you press [OK] on [Information: Do you want to change the information of the patient?], the information of the currently registered patient is shown in the patient information tab on top of the screen. Press the [X] button to maintain the current information.

NOTE

You can enter patient information using a barcode reader. Place the cursor on each item in the patient information screen and scan the barcode. Information is entered automatically.

Basically, you can use any kind of barcode reader. However, since the default setting of the input method of each product may be different, check the method supported by Bionet before using them.

- Entry methods supported by Bionet products: International standard and USB
- The products that have been tested by being connected to Bionet equipment are listed below.

No.	Manufacturer	Product Name
1	Data Logic	DL6000
2	ZEBEX	Z-3190
3	Techscan	TSK-2000

CAUTION

Each barcode reader has a product-specific initialization code. Be sure to read the user manual of the product and check the entry method before initializing it.

2) Editing a Patient

Change the information of the currently set patient.

- ① Touch [Patient] on the menu.
- ② Select [Edit].
- ③ Edit a patient information except the ID.
- ④ When you click the item, a text box where you can enter text. Press [OK] to register the items.
- ⑤ Information on the currently registered patient is shown in the patient information tab on top of the screen.

NOTE
Patient list saves up to 99 patients that have been entered so far. If the number of patients exceeds 99, the oldest patient in the list is removed.

NOTE
The patient information which has been entered or modified in the monitor does not synchronize with those in the FC central.

Chapter 5. Measuring FHR Using US

FHR is measured by obtaining the heartbeat sounds of the fetus using ultrasonic Doppler effects and subsequently calculating and saving the heart rate per minute in real-time. Since ultrasounds are reduced considerably in the air, apply sufficient amounts of ultrasonic gel to the surface of US probe to remove any air layer before using it.

FC1400 fetal monitoring device adjusts the volume automatically according to the amount of input data. When no input is detected while probe is connected, the noise will increase. When not measuring while FC1400 is running, mute it by pressing and holding the volume down key. Then when measuring, press and hold the volume up key to return back to the previous setting for the measurement.

1) Measuring FHR

■ Connecting the US probe

Connect the US probe to the US1 and US2 connection terminals on the left side of the monitor.

■ US probe



<US PROBE>

■ Basic operations following a US probe connection

When the US Probe is not connected to the monitor, nothing is displayed in FHR area, while " - " is displayed if only the US probe is connected. When you connect the US probe to the US1 connection terminal and place it to heart position of the fetus, a numerical value is displayed in

the FHR display area.


■ How to measure FHR

- ① Place a fixing belt below the waist of the patient.
- ② To remove air bubbles between the abdomen of the patient and the surface of the US probe, apply a sufficient amount of ultrasonic gel to the US probe.

NOTE

Care should be taken not to scan over a wound or incision to avoid contamination and infection.

- ③ Feel the abdomen of the patient to find the back side of the fetus and place the US probe there. If the fetus faces laterally, place the US probe on the position shown in the figure below.



NOTE

If the US Probe is placed on the side of the chest of the fetus instead of the back side, ultrasounds will not be accurately shot to the heart of the fetus; thus, the heartbeat sounds of the fetus may be frequently missed.

- ④ Move the US probe little by little around the patient's abdomen to place the US probe in a position where: the sounds from the heart of the fetus are heard relatively loud and clear;


the heart shape in the FHR display area turns green; and it blinks according to the heartbeats of the fetus. adjust the speaker volume so that the sounds from the heart of the fetal becomes appropriately louder.

- ⑤ Put the button on the upper part of the US probe into the hole in the fixing belt to fix it to prevent slipping down.
- ⑥ Put the button on the upper part of the US probe into the hole in the fixing belt to fix it so that it does not slip.
- ⑦ It takes around 2~8 seconds for FHRs to be calculated and displayed.

NOTE
If you fix the probe cable toward the patient's head, you can prevent damage to the cable and also prevent the probe from moving or shaking.

2) US screen

The US area on the screen consists of the following:



- ① Probe No.
- ② FHR Volume Icon: The current state of the Volume
- ③ Alarm range: The alarm ranges of parameters
- ④ Beat indicator: When FHR is detected as beats, the heart icon blinks.
- ⑤ FHR value: FHR indicator

3) Setup

If you select and touch the US area on the screen, the following selection window pops up. Touch the US1 area or US2 area to setup US1-related items and US2-related items, respectively.


US1		X	US2		X		
Volume	5	Alarm	On	Volume	5	Alarm	Off
		Alarm Limit		Offset	+20	Alarm Limit	
		Alarm Level	High			Alarm Level	High
		Alarm Delay	20Sec			Alarm Delay	20Sec
Signal Loss Alarm	On	Color	Yellow	Signal Loss Alarm	On	Color	Aqua

■ Adjusting sound volume

Touch the [Sound Volume] menu. Select the value to set from Off or at levels 1~9.

You can use the control panel to adjust the FHR sounds. Touch the [Sound Volume] button and select a volume from the popup window.

Also, you can directly go to the Sound Volume window using the speaker icon in the US area.



■ Alarm

You can adjust the FHR alarms for US1 and US2 separately.


Set whether to use FHR alarms.

■ Alarm Level

Set the alarm levels for FHR. You can set the alarms for patient condition alarms to Medium and Low. Refer to the alarm section for details.

■ Alarm limit

Touch [Alarm Limit] within FHR setting menu to set the alarm range of FHR. Enter numbers or touch the arrow keys to set the limit.



■ Alarm delay

Set this function to prevent alarms from ringing when FHR values exceed the alarm range but only for a short time. To set up an alarm delay, touch the [Alarm Delay] in the FHR setting menu. To not to use this function, select Off.

■ Color

Set the FHR values and Trace line colors.

■ Offset

To make it easy to distinguish when the twin fetal FHR values are similar to each other, in Trace mode and printing, the US2 value is displayed at a position where the offset value is added.

■ Signal loss alarm

Set whether to occur an alarm when the fetal heart rate is not detected. Refer to the alarm section for details.

Chapter 6. Measuring the Uterine Contraction Externally

Use external attachment-type pressure sensors to measure UA. When you attach a TOCO probe to the abdomen of the patient, it measures the relative pressures that vary with the uterine contraction of the patient, thereby recording the uterine activity.

1) Measuring Uterine Contraction

■ Connecting the TOCO probe

Connect a TOCO probe to the TOCO connection terminal on the left side of the monitor.

■ TOCO Probe




<NEW PROBE>

■ Basic operations with TOCO probe connection

When there is no TOCO probe connected to the monitor, "---" is displayed in the UC display. When you connect a TOCO probe to the monitor, a numerical value is displayed, indicating that it is ready to measure uterine activities. When you view FMs (Fetal Movements), "TOCO+FM" is displayed.

■ Measuring UA

- ① Place a fixing belt below the back of the patient.
- ② Place the TOCO probe on the highest peak of the abdomen of the patient (Fundus: located around 10cm above the umbilicus) or the place that hardens first in the abdomen of the patient.
- ③ Put the fixing button protruding on top of the TOCO probe into the hole of the fixing belt to fix the TOCO probe. At this time, the degree of tightening with the belt should be between 20 and 90 TOCO value.
- ④ Perform Zero calibration to set up a reference value.



NOTE


If a TOCO probe is connected to the monitor but not used, unreliable values may appear in the TOCO display area.

NOTE

If the type of probe (Gray or White) is not identical to the value set in FC1400, the TOCO value is not be measured. To change the probe type, check the TOCO New probe setting first.

2) UC screen

The TOCO area on the screen consists of the following:



① Title

TOCO: Displays only uterine contraction but not fetal movements.

TOCO+ FM: Displays both fetal movements and uterine contraction at the same time.

② UC Zero value: Displays the reference value that appears when Zeroing has been done.

③ UC: The degrees of uterine contraction

3) Setup

If you select and touch the UC area on the screen, the following selection window pops up.

TOCO		X
Zeroing		
Offset	10	
Auto Zeroing	Off	
Fetal Movement	FMD	
	Color	Fuchsia

■ Zeroing

Adjust the zero point of UC. If you touch it, the UC value at the moment is set as Offset value.

■ Offset


Set a reference value to use when adjusting the zero point of UC among 0, 10, and 20.

■ Auto zeroing


Initiate automatic adjustment of the zero point when a UC value is maintained below 0 for longer than 10 seconds.

■ Fetal movement

Set whether to output fetal movement-related graphs to the screen or through the printer. If the value is FMD+FM or FMD, the title of UC is shown as TOCO+FM (if Off, TOCO).



Fetal movements are automatically detected and shown as Fetal Movement graph, and as FMDs in the form of points.



■ UC Color

Set the colors of UC values and Trace.

Chapter 7. Measuring Fetal Movement

Fetal Movements can be measured in two methods: one that is automatically detected from the data input by ultrasound and one that is measured when the patient presses the marker jack when she feels the fetal movement.

■ How to use marker jack

The marker jack records the moment of fetal movement by the patient pressing it when she feels it. When the patient presses the marker jack, "|" is displayed in the FHR waveform display area. It is printed as ↑ on the recording paper.

■ Marker jack




■ Automatic measuring of fetal movements


The automatic fetal movement measurements extract information in proportion to the sizes and durations of fetal movement from the received ultrasonic Doppler signals and display the data along with uterine contraction information. In addition, if any of the obtained values is larger than the size threshold of fetal movements that has been set, the values are marked as points in the UC waveform display area. Refer to **Chapter 6. Measuring the Uterine Contraction Externally** on how to set up the threshold.

■ Setting the marker jack sound

To turn On/Off the sound of pressing Marker jack, touch the Main Menu -> System -> Marker Sound.




< LCD Display >



< Print out >

Chapter 8. Connecting the Stimulator

When you connect the stimulator to the connector and press it, a stimulator symbol appears on the screen.



WARNING	
	For a stimulator, you must use either GE products or GE-compatible products whose stability have been confirmed.

Chapter 9. Clinical Marks and Notes


To record clinical mark or notes in real-time, use [Clinical Mark] key.

Press the [Clinical Mark] key shortly to display the mark on screen or include it in Trace range printing or real-time printing. The Clinical mark is displayed in Orange triangle on screen as well as in print out.

When you press the [Clinical Mark] key long, Note window pops up. Select one of the lists to include the related information in Trace printing and real-time printing.




: Clinical Mark Key



■ Entering notes

When you press the [Clinical Mark] key long, Note window pops up. Touch the applicable item from the list and press [OK] to output the Note text in Trace printing and real-time printing.



■ Adding notes

Note lists can be added and edited only by users having Administrator rights. Touch the Main Menu -> System -> Edit Note, enter the admin password, and then use the New/Modify/Delete buttons to set the note to suit your environment.

You can only save up to 100 additional Notes, with up to 40 characters in one Note. If you try to add more than 100 Notes, an error message saying that it cannot be entered is displayed.

NOTE


The preset 15 lists cannot be edited nor deleted.

For preset password, refer to **Chapter 12, Alarm and Preset**.

Chapter 10. Trace

View the saved measurement results using the Trace function. Data measured for 72 hours are saved.

1) Trace Area



- ① This is data in the range of 30 to 240 bpm, and displays the Trace of the heart rate of the fetus.
- ② FMD: Automatically detected fetal movements are displayed.
- ③ TOCO of 0~99 units is displayed. Also, a graph of fetal movement detected by the US is displayed.
- ④ Scroll bar: Use it to move the Trace point. To navigate, touch the scroll bar or use the volume control key.
- ⑤ UC zeroing: It indicates that UC zeroing is done. The number above is the zeroing number.

⑥ Grid: Grid criteria

2) Trace

■ Trace


To view Trace data, touch [Trace] on the main menu or press the [Trace] key.



: Trace Key


■ Changing the trace time

To change the time of the Trace data, move the scroll bar of the Trace, touch the left or right side of the Trace data on the LCD, or use the volume control keys. At the bottom of the Trace screen, date information and patient name of the first Trace are displayed. If there is a Start marker within 30 seconds, the next patient's information is displayed.



By using the Volume UP key, if you press the US1 Volume UP key, the view is moved to 1 min earlier; if you press the US2 Volume UP key, the view is moved to 1 min behind.

Also, by using the Volume Down key, if you press the US1 Volume Down key, the view is moved to 6 sec earlier; if you press the US2 Volume Down key, the view is moved to 6 sec behind.



You can drag or touch the LCD screen to move as shown in the arrow above. It moves by one

page (approximately 10 minutes) when touched. ■ Exiting the trace screen

To view real-time data in Trace screen, touch the [Main Home] in the Trace menu or the [Screen Switching and Function] key on the control panel.



: Screen Switching and Function key


3) Trace Menu

Touch the menu in the Trace screen. The Trace menu appears.

EXIT				Data Save	Print Start	Main Home

■ Saving data

Save a certain period of data from the Trace to the USB storage device. If you insert a USB storage device into the slot on the back of the monitor, a USB device recognition icon is shown on top of the screen. Once a USB storage device is recognized, the [Data Save] button is activated when you touch the main menu in the Trace.



When you touch the [Data Save] button, data is saved to USB. The first position of the Trace is saved as data image up to 3 hrs within one page. Duration depends on the print speed, 30 min for 3cm/min, 60 min for 2cm/min and 90 min for 1cm/min. If the patient is changed in the meantime, saving is suspended. If there is a Start marker within the first 30 seconds of the first position of the Trace, data is saved from that position. The following window appears, and when the progress bar reaches the end after about 10 seconds, saving is completed and the window automatically disappears.

You can find the data in the *WNewFC1400Wdata* path of the recognized storage device. The name of the save file is: *{Patient ID}_{Patient Name}_{(Gestational-age)}_{Year Month Day Hour Minute Second}.jpg*

NOTE

If the patient's ID or name contains special characters, images cannot be created.



WARNING	
	Remove the USB storage device from the slot only after data saving is completely finished and the <i>Saving data</i> window disappears. Do not remove the USB device from the monitor during the data transfer.

■ Printing

You can print a certain time of the Trace. Up to 1 hour of data is printed from the Trace position on the current screen. If the patient is changed in the meantime, data until the patient change is printed. If there is a Start marker within the first 30 seconds of the Trace, data is saved from that position.

You can also touch the [Print] key on the control panel to perform the same feature.

■ Main home

The Trace mode is terminated, and data is shown on screen in real-time.

WARNING	
	Trace data is saved periodically by every 1 minute. Thus, up to 1 minute from the end of data may not be saved. Additionally, no data is saved if you turn off the monitor within 1 minute of powering on.

WARNING



Trace data may be erased after the SW upgrade.

Chapter 11. Printing

FC1400 fetal monitoring device offers 2 kinds of printing: real-time printing and Trace printing. In the printing process, print icon is seen on the right top. There are 2 types of Print icons: Printing and Not printing.



(Printing)



(Not printing)

To start or stop print, press the following [Print] key.

To start or stop normal printing, press the key shortly. To accelerate the printing, keep pressing the key. If you take your hand off the key during the accelerated printing, the maximum of 2.5 cm is printed more and printing stops.



: Print start/stop key

1) Real-time printing

To print currently inputted data, press the [Print] key or touch the [Print Start] in the menu. If any printing is in progress, a print icon is displayed on the upper right side of the screen. To stop printing, press the [Print] key again or touch the [Print Stop] in the menu.


When printing in real-time has been stopped, paper is fed in a certain length to facilitate cutting.

Printing in real-time is done in any of three speeds: 1, 2, 3 cm/min.

■ Changing the speed of real-time printing

To change FC1400's speed of real-time printing, touch the Main Menu -> Printer. Touch the [Print Speed] in the menu and select the desired speed to change the speed.

Printing speeds is not changed while the real-time printing is in progress. To change the speed, stop the printing first, change the speed, and then resume printing.



■ Changing the type of printing paper

To change the type of paper, touch the Main Menu -> Printer. Touch [Scale] in the menu and select the desired type of paper.




■ Changing the intensity of line

To change the intensity of print line of paper, touch the Main Menu -> Printer. Touch [Print line] in the menu to change the line intensity.

■ Changing the size of printing paper

To change the size of printing paper, touch the Main Menu -> Printer. Touch [Size] in the menu and select the desired size of printing paper between B5 and A4 in the following menu.



■ Printing the grid line

To print grid line on the fax paper, touch the Main Menu -> Printer -> Grid. Set the Grid to On to print the grid on the paper.

2) Printing the Trace

To print the Trace data, move to the point from which you want to print on the Trace window and press the [Print] key or touch the [Print Start] in Trace menu. The Trace is printed up to the next measurement or up to 1 hour. The Trace printing speed is 30cm/min (when NST print speed is 2cm/min or 3cm/min) or 20cm/min (when NST print speed is 1cm/min).

3) Loading the Paper

If there is no paper during printing, *No paper* alarm occurs. In that case, load paper or press [Alarm silence] key to release the alarm. To load the paper, open the printer cover by pushing back the Printer cover open lever on the right side of the monitor, and load paper with printable side up and close the cover.

4) Registering hospital name

To register hospital name on printing paper, touch the Main Menu -> Printer -> Hospital.

Chapter 12. Alarm and Preset

Alarms are largely divided into patient condition alarms and product condition alarms (INOP Alarms).

1) Patient Condition Alarms

Patient condition alarm occurs when the value exceeds the maximum and minimum limit values for alarms. There are two levels of alarms: Medium and Low, which differ in terms of the order of ringing and volume.

Medium alarms are shown as ** on the alarm list when they occur, and Low alarms, as *.

The FHR alarm occurs when the value is out of the normal range during the set Alarm Delay time.

If more than one alarm sounds, alarm messages appear successively on the alarm condition window.

If two or more alarms occur, the alarm at the highest level will sound. If an alarm has already occurred, it will sound at the previous level even if you change its level. To immediately apply the alarm to the changed level, turn it off and on to trigger the alarm again.

MEDIUM	A low-pitched sound will be repeated three times and paused for a while.
LOW	A low-pitched sound will be repeated two times and paused for a while.

2) Product Condition Alarms

Product condition alarms are used to indicate that FC1400 is not operating properly and its capability of detecting the dangerous conditions of the patient is not reliable.

Alarms include those related to connector connection and technical alarms such as *FHR disappearing, No paper, and Battery Low*.

"Ding-dong"	When the connector is being disconnected
Technical Alarm	A low-pitched sound is repeated one time and paused for a while.

Alarm Level	Volume	Alarm Interval	Beep Rate	Sound Pressure level[dB]
Medium Level	1 step	Every 8 sec	3 beeps	46.2 ~ 53.1
Low Level	1 step	Every 8 sec	2 beeps	46.1 ~ 53.0
Medium Level	5 steps	Every 8 sec	3 beeps	63.4 ~ 68.2
Low Level	5 steps	Every 8 sec	2 beeps	63.1 ~ 67.6

NOTE

High-priority alarms indicate that immediate operator response is required. Low-priority alarms indicate that operator awareness is required.

3) Visual Alarms

Alarm messages appear in the Alarm State window on the screen. If more than one alarm occur, messages are changed every 2 seconds. The * marks of the alarm messages match the level of the alarm. Medium Alarms are marked with **, Low Alarms with *, and product condition alarms with none.

When a patient condition alarm occurs, the color of the numerical value in the numeral window is changed to red.

Visual Alarms are maintained even after you change the alarm state into Alarm Silence or Alarm Pause.

No Paper and *Battery Low* alarms are displayed in the alarm state window.

Visual alarm		
Alarm priority		LED Flashing Duty Cycle
PATIENT	MEDIUM	** Indicates Medium Alarms

CONDITION		Yellow LED; Flickers once a second
	LOW	* Indicates Low Alarms
		Yellow LED; Stay on
PRODUCT CONDITION		Green LED; Alarm occurs due to technical alarm and stays on.

4) Alarm LED

There are two types of Alarm LED around the exterior of the monitor to help checking the alarm condition from a distance.

Yellow LED: An alarm occurs due to patient condition.

Medium Alarm: Flickers once every second.

Low Alarm: Stays on

Green LED: An alarm occurs due to Technical Alarm and stays on.

5) Audible Silence and Pause

When various alarms occur during the operation, use this function to silence or pause them to verify the message.

Press the [Alarm Silence] key once to sound off (Silence) for 1 minute. If a new alarm occurs during the Silence operation, alarm sound will ring again.

Press the [Alarm Silence] key twice to make the alarm in pause state for 5 minute. In Pause state, no alarm will sound even if a new alarm is generated.

If you press the [Alarm Silence] key 3 times, sound alarms are turned off and no sound will be made even if a new alarm occurs.

Visual alarms are displayed until the situation is resolved. To mute the alarm, touch [Alarm Silence] key on the operational panel.



: Alarm Silence key


WARNING	
	When monitoring patient conditions, do not rely solely on audible signals. Reducing or turning off the volume of audible signals may endanger the patients.

6) Alarm Latching

Audible alarms continue to ring until you take an appropriate action and the situation is resolved. Once you take measures and the alarm situation is over, both audible and visual alarm disappear.

7) Alarm History

You can check the history of alarms that occurred to the patients. Touch the Main Menu -> Alarm History in the alarm message area to view the alarm history window. If you admit a new patient or change to another patient, the Alarm History is cleared.



8) Preset

Preset values are those loaded by default on system rebooting, which can only be changed by an administrator. Volume adjusting applies to current patient only, and when you restart FC1400 or use it for another patient, it operates using the preset values. However, if you turn it on within 30 seconds after abnormal shutdown, the previous setting value is maintained assuming that it is still registered with the same patient.

Preset consists of alarm volume, alarm function usage, alarm range, alarm level, alarm delay time, US1 Volume, and US2 Volume.

To set Preset, touch the Main Menu -> Preset.


As Preset is accessible only by an administrator, you need a password in order to enter the menu.

WARNING	
	The initial Preset password in Factory mode is 1234. Change the Preset password after installing FC1400 to avoid unauthorized access. To change the Preset password, touch the Main Menu->System->Change PW and enter password. Be aware, when the changed password is lost, the only way to reuse FC1400 is initializing it to factory mode. On initializing, all data will be removed. Make sure to keep the password safe.

9) Adjusting Alarm Volume

You can adjust the Audible Alarms volume to 1~5 levels or Off.

- ① Touch the Main Menu->Preset->Alarm Volume.
- ② Select the desired volume on the popup window.



- ③ Check if the current alarm volume is displayed on the alarm icon on the screen.



: Alarm Volume condition icon Off, levels 1~5



: Alarm off icon. Turn off the alarm. (Off)




: Alarm Silence



: Alarm Pause

10) Setting All Alarms ON/OFF

You can activate or deactivate alarm for each parameter. Touch the Main Menu -> Preset-> Alarm On/Off and set the condition by pressing the button in Alarm Parameter window. To apply default values, select the [Default] button and press [OK]. For default values, see **Chapter 19**.



11) Setting Ranges of All Alarms

You can set the alarm range for each parameter. Touch the Main Menu->Preset->Alarm Limit and set the alarm range for each parameter in the Alarm Limit window. To apply default values, select the [Default] button and press [OK].

The screenshot shows the 'Alarm Limits' window with a blue header and a close button (X). The window contains two columns for parameters US1 and US2. For each parameter, there are 'High' and 'Low' settings. The 'High' values are set to 160, and the 'Low' values are set to 120. At the bottom, there are 'Default' and 'OK' buttons.

	US1	US2
High	160	160
Low	120	120

11) Setting Levels of All Alarms

You can set the alarm level for each parameter. Touch the Main Menu->Preset->Alarm Level and set the alarm level for each parameter in the Alarm level window. To apply default values, select the [Default] button and press [OK].

The screenshot shows the 'Alarm Levels' window with a blue header and a close button (X). The window contains two columns for parameters US1 and US2. For each parameter, there is a level setting. The level for both US1 and US2 is set to 'Medium'. At the bottom, there are 'Default' and 'OK' buttons.

	US1	US2
Level	Medium	Medium

13) Default Setting

You can set default values in Alarm Parameter, Alarm Limit and Alarm Level. Touch the Main Menu->Preset->Default Setting to apply default values in the Preset settings. For default values, see **Chapter 19**.

NOTE
When using the monitor through FC central, all setting values are synchronized once you are out of the main menu and parameter window.

14) Signal Loss Alarm

Signal Loss Alarm occurs when FHR is lost for a certain period.

- 100 % Signal loss: Absent FHR for last 75 seconds.
- 70 % Signal loss: Less than 30% acceptable FHR for last 5 minute.
- 65 % Signal loss: Less than 35% acceptable FHR for last 10 minute.

The alarm is displayed as a Product condition alarm.

The signal loss alarm is enabled/disabled on the US1 and US2 menu.

Chapter 13. Network

1) Setting Network

To set up a network, touch the Main Menu -> Network.

There are 5 setting criteria.

EXIT	IP	Central Server	Wireless :Off	Central :Off	Bed Num :7	
Up						

- IP: Monitor network information

- Central IP: Central server information

Enter the IP of the PC where the remote PC (central) program is installed, with fixing the Port to 5000.

- Wireless setting

If using the Wireless LAN, insert the wireless LAN module into the USB slot on the rear side, and touch the Wireless menu to set it to On. Set Central to ON. When the both menu above are set to ON, the wireless LAN function is activated and network menu configuration looks different that appears when you touch the IP button. When the wireless is Off, the network is connected using wired LAN configuration. For more information about FC Central network configuration using wireless, refer to the FC Central Installation Guide.

- Setting Central

If setting Central to On, you can connect the monitor to the Central. If you set it to Off, the connection is terminated.

WARNING	
	Before setting the Wireless and Central to ON, check the Wireless LAN module is properly installed in the monitor. Wireless is not selectable unless it is well recognized.

WARNING	
	VLAN Network If data is exchanged within a single network, an independent VLAN network for the clinical information system must be established. A network system that detects and defends against denial-of-service attacks must be established through the installation of equipment dedicated to DDos defense. When using wireless, ensure a proper AP security; if it not available, a wired connection is recommended.

Mandatory Action	
	DO NOT remove the USB storage device until it is recognized completely.


- Changing the bed number

To change the bed number, which is required for monitor identification, touch Bed No. in the setting menu. You can select the number from 1~16

NOTE
Wireless AP's ID are recognized only if it consists of alphanumeric characters.

Chapter 14. General Settings

Go to System in the main menu to define the general settings.



The menu of the system is as follows:

EXIT	Date	Time	Lang :English	Version	Touch Tone :3	Demo :Off
Up	Edit Note	Marker Sound:On	Factory	Change Admin PW	Volume range: 3	Protocol Version:1.0

1) Changing the Date

To set up a date, touch System --> Date.

2) Changing the Time

To set up time, touch System --> Time.

3) Changing the Language

To set up language, touch System --> Language.

4) Checking the System Version

To check the version of the system, touch System --> Version.

FC1400 Device Information		X
S/W Version		
UI		1.0.0
F/W #1		1.0.0
F/W #2		1.0.0
H/W Version		
Digital		0
Analog		0
H/W Address		00:00:00:00:00:00

S / W Version: Software Version

UI : Version of UI software

F/W #1 : Version of fetal module firmware

F/W #2 : Version of printer module firmware

H / W Version: Hardware Version

Digital : Digital board version

Analog : Analog board version

H / W Address : Mac Address

5) Changing the Touch Tone

To change the volume of touch sound, touch System --> Touch Tone.

Set Off or select 1~5.

6) Demo Operation

To operate FC1400 in the demo mode, touch System --> Demo and set it to On.

7) Editing Note

Touch it for New/Edit/Delete of Note list. This function can be only used by an administrator. When correct Preset password is inserted, setting window pops up. For more information, see **Chapter 9, Clinical Marks and Notes**.

8) Marker Sound

Touch it to turn On/Off the sound when clicking a Marker. For more information, see **Chapter 7, Measuring Fetal Movement**.

9) Factory Mode

Only the administrator can enter the Factory mode. In order to change the item, request for a service to Bionet Headquarter or its representative.

10) Changing Admin PW

Change Admin password in this menu.

At first, enter the current Admin password correctly. Once it is confirmed, you can proceed to change it by entering a new and confirm password.

The default Admin password is **1234**.

WARNING	
	Admin password cannot be recovered after a change. If password is lost, reset to Factory mode is the only way to reuse FC1400, and all the measured data and set values will be lost. Make sure to keep the password safe.

11) Volume Range

Use this menu to change volume range. 7 is the loudest and 0 is the quietest volume.

12) Screen Output Mode

The screen supports two output modes: Graphic Mode and Text mode. To change a screen output mode, touch the [screen switching and function] key on the control panel in the main menu.



: Home and function key

13) Protocol Version

Set the version of protocol according to the version of FC Central. To use Note function in FC1400 and have it synced to FC Central, select protocol version 1.2 and FC central version 1.2.2 or higher.

Chapter 15. Non Stress Test (NST)

The Non-Stress Test Timer measures the NST time that passed. You can print the only the data measured during a certain period of time.

1) Measuring NST

■ Starting NST

Press the [Print] key to start NST.

When the set NST time has passed, alarms will sound to indicate the NST is over.

To mute the NST ending alarm, touch the [Alarm Silence] key on the control panel.

■ Ending NST

Press the [Print] key again to stop the NST even if the NST measurement time is not over. The printing stop as well.

2) Setting NST

Touch the Main Menu -> Setting -> Printer-> NST causes the following setting window to pop up.

NST Time			
Off	10Min	20Min	30Min
40Min	50Min	60Min	90Min

■ NST Duration

Set the NST duration. Select OFF, 10, 20, 30, 40, 50, 60, and 90 min.

Chapter 16. CTG

CTG is FC1400 fetal monitoring device' interpretation about the Cardio-TOCO Gram. The analysis begins with the data collected from that moment when you press the [Printer] button.

1) Setup

▣ Setting CTG

Touch the Main Menu → Printer → CTG. Set the CTG Print to On to print out the CTG result.

▣ Turning off CTG

Touch the Main Menu → Printer → CTG. Set the CTG Print to Off.

2) Printing CTG

▣ Starting CTG


CTG analysis starts when you start the real-time printing. The CTG analysis result requires equal to or more than 10 minutes printing.

▣ Stopping CTG

When you stop the real-time printing, the CTG final result is printed after 5 seconds. To print out CTG final result, you have to measure data at least 10 minutes.

3) CTG Measurement results

After CTG is measured, the result is printed out as below.



4) Glossary of CTG Terms

■ Base FHR / Average BaseFHR

Excluding constant or irregular FHR, significantly changed FHR, or sections where the difference from the baseline exceeds 25 bpm, it must run over 2 minutes in increments of 5 bpm with an average FHR value over 10 minutes.

■ Number of UC

It is number of uterus contractions during measurement.

■ High Episode / Low Episode

If the changing width of the last 5 minutes is less than 30msec, it is a Low Episode, and if it is more than 32ms, it is a High Episode.

■ Short Term Variability

Short Term Variability (STV) is the beat-to-beat differences between consecutive heart beats. If the value of STV is 2.6 ms or less, carefully monitor the condition of the fetus and perform additional tests if the fetus is judged to be in critical condition.

■ Acceleration


It is visually apparent increase (onset to peak is < 30sec) of FHR above baseline. Peak is ≥ 15 bpm. Duration is ≥ 15 sec and <2 min.

■ Smaller Acceleration

It is visually apparent increase (onset to peak is < 30sec) of FHR above baseline. Peak is ≥ 10 bpm and <15 bpm. Duration is ≥ 15 sec and <2 min.


■ Late Deceleration

It is visually apparent gradual decrease (onset to nadir is ≥ 30 sec.) of FHR below baseline. Return to baseline is associated with a uterine contraction. Nadir of deceleration occurs after the peak of the contraction. Generally, the onset, nadir and recovery of the deceleration occur after same time as the onset, peak, and recovery of the contraction.




■ Early Deceleration

It is visually apparent gradual decrease (onset to nadir is ≥ 30 sec.) of FHR below baseline. Return to baseline is associated with a uterine contraction. Nadir of deceleration occurs at the same time as the peak of the contraction. Generally, the onset, nadir, and recovery of the deceleration occur at the same time as the onset, peak and recovery of the contraction.



■ Variable Deceleration

It is visually apparent abrupt decrease (onset to nadir is $<30\text{sec}$) in FHR below baseline. Decrease is $\geq 15\text{ bpm}$. Duration is $\geq 2\text{ min}$ and $< 10\text{ min}$.



■ Fetal Movement per Hour

Fetal movement - Indicates the fetal movement as a number of times per hour.

Chapter 17. Message List

1) Patient Alarms

Short Msg	Figure Msg	From
FHR* Medium	FHR* xxx < yyy	FHR*
FHR* LOW	FHR* xxx>yyy	FHR*

2) Technical/INOP Alarms

Short Msg	Figure Msg
No Paper	Main
Low Battery	Main
Signal Loss	FHR* 100% Signal Loss FHR* 70% Signal Loss FHR* 65% Signal Loss

Chapter 18. Simple Troubleshooting

1) Troubleshooting and Solutions

- ① If the probe comes off during operation, a "---" sign is displayed, and a "Ding" sound is heard. In that case, check the probe's connection condition and connect it again.
- ② If paper has run out during operation, a *No Paper* sign is displayed on the LCD screen. In this case, open the printer, check if the recording sheets have run out, replenish the record sheets, and close the printer.

WARNING	
	If Touch is not calibrated properly, FC1400 may not operate properly. You must calibrate the touch input as written in the operation manual.

2) Performing Periodic Inspections

Just like all kinds of medical equipment, conduct safety inspections on FC1400 periodically (once a year). Refer to the service manual provided by Bionet for inspection items.

3) Cyber Security Issues

1. If equipment is stolen or lost, immediately report it to the hospital staff or manufacturer. Upon receipt of a report, the hospital network administrator must take measures to prevent the equipment from accessing the hospital network.

2. If a cyber security threat is detected while using the equipment, immediately disconnect it from the network and contact the hospital staff or manufacturer.

※ For manufacturer contact information, refer to the table of contents of How to Contact Us.

Chapter 19. Product Specifications

General specification

Dimension	296(w)x.305.5(H)x97.5(D)(Approx. 2.9Kg)
Display	7" Wide (800 X 480)
Recorder	Method : Thermal Array Print Type : Roll type Print speed : 1,2,3cm/min,(Real time) 30cm/min (Trace, 2,4 cm/min setting) 20cm/min (Trace, 1 cm/min setting) Paper feeding function

Performance Specification

Fetal Heart Rate	Input Signal : Ultrasound pulsed doppler FHR detection method : Autocorrelation FHR range : 50~210 FHR accuracy : 120~160 : ± 1 bpm Except 120~160 : ± 2 bpm
Ultrasound Transducer	Operating mode : PWD Mode Transducer Type : 7-crystal Ultrasound frequency : 1.0MHz Pulse Repetition Frequency : 3125Hz Spatial-Peak Temporal Average Intensity : $<10\text{mW/cm}^2$
Uterine Contraction	Input Source : External Transducer Reference Control : One touch switch Auto zeroing Measurement range : 0~99
Auto CTG Analysis	Average Baseline FHR Number of TOCO Number of Acceleration Number of Deceleration : Late, Early, Variable High/Low Episode Short term Variability Signal Loss * CTG Analysis results is printed out every 10 minutes
Data Storage	Storage for 72 Hours

Power Specification

Power	DC Input: 18V $\overline{\text{---}}$, 2.8A / Adaptor "USE THE ONLY BridgePower Corp BPM050S18F02"
Battery(Optional)	Li-ion : 4 hours (charging) 2 hours(discharging)
External Link	LAN, Wi-fi, USB,SD

Default Alarm Setting

Alarm Parameter	US1/ US2	ON	
Alarm Limit	US1/US2	160	120
Alarm Level	US1/US2	Medium	

■ Intended Use

FC1400 detects and displays twin Fetal Heart Rate (FHR), Fetal Movement (FM), and Uterine Activity (UA) in real-time on the color LCD viewer, and also provides the fetal heart beat sound with internal speaker. 12 hours of tracing may be saved and later retrieved for printing. It intended to be used by trained healthcare personnel. It is not intended for home use.

■ Indication for Use

FC1400 is suitable for use when there is a need to monitor the following physiological applications:

- Single, twin or triplet fetal heart rates by means of ultrasound
- Uterine activity - externally sensed
- Fetal movement - maternally sensed and externally via ultrasound

■ Contraindications:

Do not use FC1400 to directly monitor patients during water births, in whirlpool or submersion water baths, during showers, or in any other situation where the mother is immersed in water.

Doing so may result in electrical shock hazard.

Appendix A.

Manufacturer's declaration - electromagnetic immunity

The FC1400 system is intended for use in the electromagnetic environment specified below. The customer or the user of the FC1400 system should assure that it is used in such an environment			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic Environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV Contact 8 kV Air	6 kV Contact 8 kV Air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast Transient / burst IEC 61000-4-4	2kV for power supply lines 1kV for input/output lines	2kV for power supply lines 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	1 kV differential mode 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) Magnetic field IEC 61000-4-8	3.0 A/m	3.0 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Voltage dips, short Interruptions and Voltage variations on power supply input lines IEC 61000-4-11	<5% U_t (>95% dip in U_t) for 0.5cycle 40% U_t (60% dip in U_t) for 5 cycle 70% U_t (30% dip in U_t) for 25 cycle <5% U_t (<95% dip in U_t) for 5 s	<5% U_t (>95% dip in U_t) for 0.5cycle 40% U_t (60% dip in U_t) for 5 cycle 70% U_t (30% dip in U_t) for 25 cycle <5% U_t (<95% dip in U_t) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BM1 system requires continued operation during power mains interruptions, it is recommended that the BM1 system be powered from an uninterruptible power supply or a battery
Note: U_t is the a.c. mains voltage prior to application of the test level.			

The **FC1400** system is intended for use in the electromagnetic environment specified below.
The customer or the user of the **FC1400** system should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the BM1 system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$

Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	Recommended separation distance $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (a) Should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1) U_t is the A.C. mains voltage prior to application of the test level. Note 2) At 80 MHz and 800 MHz, the higher frequency range applies. Note 3) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the EUT is used exceeds the applicable RF compliance level above, the EUT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the EUT.			
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V / m.			
Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the FC1400 system.			
The FC1400 system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the FC1400 system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FC1400 system as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power (W) of transmitter	Separation distance (m) according to frequency of transmitter		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.70	3.70	7.37
100	11.70	11.70	23.30
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

Immunity and Compliance Level			
Immunity test	IEC 60601 Test Level	Actual Immunity Level	Compliance Level
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms, 150 kHz to 80 MHz	3 Vrms, 150 kHz to 80 MHz
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m, 80 MHz to 2.5 GHz	3 V/m, 80 MHz to 2.5 GHz

Guidance and manufacturer's declaration - electromagnetic immunity

The **FC1400** system is intended for use in the electromagnetic environment specified below.

The customer or the user of the **FC1400** system should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	FC1400 system must be used only in a shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location with a minimum RF shielding effectiveness and, for each cable that enters the shielded location
Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Note 2) It is essential that the actual shielding effectiveness and filter attenuation of the shielded location be verified to assure that they meet the minimum specification.

a- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength outside the shielded location in which the EUT is used exceeds 3V/m, the EUT should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as relocating the EUT or using a shielded location with a higher RF shielding effectiveness and filter attenuation.

Appendix B. Maintenance, Care, and Service

1) Mechanical hazard

WARNING	
	Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use with care when handling and protect from damage when not in use. DO NOT use a damaged or defective probe. DO NOT drop the probes or subject them to other types of mechanical shock or impact.

WARNING	
	A defective probe or excessive force can cause patient injury or probe damage: • Observe depth markings and do not apply excessive force when inserting or manipulating intercavitary probes. • Inspect probes for sharp edges or rough surfaces that could injure sensitive tissue. • DO NOT apply excessive force to the probe connector when inserting into the probe port. The pin of a probe connector may bend.

2) Biological hazard

WARNING	
	To avoid the risk of disease transmission: • Must use protective barriers (gloves and probe sheaths). Follow sterile procedures when appropriate. • Thoroughly clean probes and reusable accessories after each patient examination and disinfect or sterilize as needed. • Follow all infection control policies established by your office, department or institution as they apply to personnel and equipment.

3) Electrical hazard

WARNING	
	In case Gel contacts internal electronic device, defective probe may cause electrical shock. Prior to each use, visually inspect the probe lens and case area for cracks, cuts, tears, and other signs of physical damage. DO NOT use a probe which appears to be damaged until you verify functional and safe performance. Perform a more thorough inspection, including the cable, and connector, each time you clean the probe. DO NOT kink, tightly coil, or apply excessive force on the probe cable. Insulation failure may result. Warning: To avoid risk of electric shock, this equipment must only be connected to a supply main with ground Do not modify this equipment without authorization of the manufacturer If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment. Do not touch signal input, signal output or other connectors, and the patient simultaneously. Refer servicing to qualified personnel of Bionet Co., Ltd. Power supply is specified as a part of ME Equipment.

4) Probe acoustic output hazard

WARNING	
	Ultrasound can produce harmful effects in tissue and potentially result in patient injury. Always minimize exposure time and keep ultrasound levels low when there is no medical benefit.

5) Probe head: Waterproof

MANDATORY ACTION	
	From US probe bottom to 2~3 cm, waterproof IPX is possible. Do NOT immerse the probe bottom into any liquid beyond 2~3 cm from probe bottom. Never immerse the probe connector into any liquid.

Appendix C. Ultrasound Power

Use of Diagnostic Ultrasound

The American Institute of Ultrasound in Medicine (AIUM) has published a document entitled "Medical Ultrasound Safety".

This three-part document covers Bioeffects and Biophysics, prudent Use and Implementing ALARA.

Ultrasound users should read the AIUM documents to become more familiar with Ultrasound safety. A copy of this document is included as part of the documentation package (Document 2163920-100).

AIUM

14750 Sweitzer Lane

Suite 100

Laurel, MD, USA 20707-5906

telephone 1-800-638-5352.

In accordance with US FDA Guidelines, the overall maximum acoustic SPTA intensity for the product is limited to 100 mW/cm² and MI is limited to 1.0.

1) Measurement Precision and Uncertainty

	Center Frequency	Acoustic Power	Peak Rarefractional Pressure	Acoustic Intensity
Measurement Uncertainty	± 2 %	±5%	±15%	± 25%

2) Maximum Output Summary

Operating Mode	DOP Probe
Pulsed Doppler Mode	O

3) Maximum Probe Temperature (Degrees C)

Probe	Max Temperature		
	With TMM Phantom	In Air	Mode
US	33.2	21.6	PWD Mode

Lens temperature monitored for 30 min.

Measurement uncertainty: +/-0.5-degree C.

Ambient temperature: 22.1 degree C

4) Table Key

IEC	FDA	Meaning IEC60601-2-37 / FDA&NEMA UD2,UD3
α	a	Acoustic Attenuation Coefficient / Derating factor (usually 0.3 dB/cm-MHz)
Aaprt	Aaprt	-12db Output Beam Area / Active aperture area
CMI	-	Normalizing Coefficient
Deq	Deq	Equivalent Aperture Diameter
d-6	d-6	Pulse Beam Width / Beam diameter at -6 dB
deq	deq	Equivalent Beam Diameter
f_{awf}	fc	Acoustic Working Frequency / Center frequency
Ipa	Ipa	Pulse-Average Intensity
Ipa, α	Ipa.3	Attenuated Pulse-Average Intensity
Ipi	PII	Pulse-Intensity Integral
Ipi, α	PII.3	Attenuated Pulse-Intensity Integral
I _{ta} (z)	ITA	Temporal-Average Intensity
I _{ta} , α (z)	ITA.3(Z)	Attenuated Temporal-Average Intensity at depth z
I _{zpta} (z)	ISPTA(Z)	Spatial-Peak Temporal-Average Intensity
I _{zpta} , α (z)	ISPTA.3(Z)	Attenuated Spatial-Peak Temporal-Average Intensity
MI	MI	Mechanical Index
P	Wo	Output Power / Time average acoustic power at the source
P α	W.3(Z)	Attenuated Output Power / Time average acoustic power derated to depth z
P1	Wo1	Bounded Output Power / Power emitted from the central 1cm of aperture
pi	PII	Pulse Pressure Squared Integral / Pulse intensity integral
pr	pr	Peak-Rarefactional Acoustic Pressure
pr α	pr.3	Attenuated Peak-Rarefactional Acoustic Pressure
prr	PRF	Pulse Repetition Rate / Pulse repetition frequency
TI	TI	Thermal Index
TIB	TIB	Bone Thermal Index
TIC	TIC	Cranial-Bone Thermal Index
TIS	TIS	Soft-Tissue Thermal Index

td	PD	Pulse Duration
X, Y	x-12,y-12	-12 dB Output Beam Dimensions
Z	Z	Distance from the Source to a Specified Point
Zb	Zsp	Depth for TIB / Depth at which the relevant index is maximum
Zbp	Zbp	Break-Point Depth
Zs	Zsp	Depth for TIS / Depth at which the relevant index is maximum

Acoustic Output Tables

MC65R1S – Pulsed Doppler Mode

Index				MI	TIS			TIB	TIC
					scan	Non-scan		Non-scan	
						Aaprt<= 1	Aaprt > 1		
Global Maximum : Index Value				0.0164842	-	0.00168143	-	0.0130577	0.00869565
	IEC	FDA	Unit						
Associated	pra	pr.3	(MPa)	0.0164807					
Acoustic	P	Wo	(mW)		-	0.4		0.4	0.4
Parameter	min of [Pα(zs), I _{ta,α} (zs)]	min of [(W.3(Z1), I _{TA.3} (z1)]	(mW)				-		
	Zs	z1	(cm)				-		
	zbp	zbp	(cm)				-		
	zb	zsp	(cm)	1.8				1.8	
	z at max. I _{pi,α}	zsp	(cm)						
	deq(zb)	deq(zsp)	(cm)					1.0865	
	fawf	fc	(MHz)	0.999572	-	0.999572	-	0.999572	0.999572
	Dim of	X	(cm)		-	0.4	-	0.4	0.4
	Aaprt	Y	(cm)		-	0.4	-	0.4	0.4
Other	td	PD	(μsec)	59.9647					
Information	prr	PRF	(Hz)	3906					
	pr at max. I _{pi}	pr@P _{II} max	(MPa)	0.0175374					
	deq at max. I _{pi}	deq@P _{II} max	(cm)					1.0865	

	Focal Length	FLX	(cm)		-	2	-		2
		FLY	(cm)		-	2	-		2
	Ipa,α at max. MI	IPA.3@MI _{max}	(W/cm ²)	0.00154839					
Operating Control Conditions	Frequency		(MHz)	1.0	-	1.0	-	1.0	1.0

Appendix D. Abbreviation and Symbol

Abbreviation and Symbol of manual or system operation is arranged in alphabetical order.

Abbreviations

		A
AC	alternating current	
		B
		C
C	Celsius	
cm, CM	centimeter	
		D
DC	direct current	
		E
EMC	electromagnetic compatibility	
EMI	electromagnetic interference	
		F
F	Fahrenheit	
		G
g	gram	
		H
HR	heart rate, hour	
Hz	hertz	
		I
Inc	incorporated	
		J
		K
kg, KG	kilogram	
		L
L	liter, left	

lbs, LBS pounds

LCD liquid crystal display

LED light emitting diode

M

M mean, minute

m meter

MIN, min minute

MM, mm millimeters

MM/S millimeters per second

MMHG, mmHg millimeters of mercury

mV millivolt

N

O

P

		Q
		R
		S
sec	second	T
Temp, TEMP	temperature	U
		V
V	volt	W
		X
X	multiplier when used with a number (2X)	Z

Symbols

&	and
°	degree(s)
>	greater than
<	less than
–	minus
#	number
%	percent
± , +/-	plus, or minus

Product Warranty

Product Name	Fetal monitoring device
Model Name	FC1400
Product Permit No.	
Date of Issuing Product Permit	
Lot Number	
Warranty Period	One year from the date of purchase
Date of Purchase	Day Month Year
Customer Information	Name of Hospital: Address: Name: Phone:
Name of Seller	
Name of Manufacturer	

- ※ Thank you for purchasing FC1400.
- ※ This product is a medical device.
- ※ This product has passed thorough quality control and strict inspection.
- ※ Compensation criteria for repair, exchange, and refund of this product are in accordance with the Fair Trade Commission Notice "Consumer Damage Compensation Regulations."

■ Customer Support and Repair Station at Head office

5th fl. 61, 31-gil Digital, Guro-gu, Seoul
(08375)

Tel : +82-2-6292-6410 / Fax : +82-2-6499-7789

E-mail : service@ebionet.com

Bionet Co., Ltd.

Model Name: FC1400